

EFFIVET SKIN & WOUND can be used to moisten, clean and debride wounds without harming healthy tissue. This ready to use spray is pH neutral and wont sting or irritate wounds. It can be used to provide relief from scratching, for hot spots and can be used to clean the ear and eye lid areas.

How to use EFFIVET SKIN & WOUND Hydrogel

- Remove excess hair from wound area.
- Apply liberally on the wound.
- Allow EFFIVET to saturate the area.
- If a dressing is required, apply EFFIVET liberally to the dressing. No rinsing is required.
- The skin may “pinken” after application due to increased blood flow to the wound site, which will speed up healing.

Precautions for use:

For external use only. Discontinue use if irritation develops and contact your veterinarian. For animal use only. Keep out of reach of children. Always consult with a veterinarian for systemic infections.

Can be licked by the animal.

Manufactured by:

Cf Pharma Ireland Ltd,
73B Hebron Industrial Estate
Kilkenny, Ireland

Marketed and Distributed by:

Ascendis Animal Health (Pty) Ltd
Unit 45B, 45 Parkview Street,
Highway Business Park,
Centurion,
South Africa

EFFIVET

SKIN & WOUND

HYDROGEL

CLEANS
PROTECTS
PREVENTS INFECTIONS

Hydrogel increases contact time,
therefore increases results.

Alcohol free
pH Neutral

Cleans and soothes
Eliminates odour
Relief from scratching
Wounds & Cuts
Ears & Eyes
Hot Spots
Insect Bites

Components:

Each bottle of EFFIVET contains:
0,024 % Hypochlorous Acid
Excipients: Ionised water, sodium chloride, sodium magnesium fluorosilicate gel and monosodium phosphate.

Did you know?

EFFIVET:

- Helps clear up hotspots
- Relieves scratching
- Can be licked or ingested by animals
- Non-irritating
- Cleans and soothes
- Doesn't sting or burn

Pack size 250 ml

Storage:

Store at room temperature, at or below 25 °C away from direct sunlight and heat. DO NOT ALLOW TO FREEZE. Do not transfer into another container. Recap after each application. Keep out of reach of children.

Bar code

See base of container for batch number and use by date.